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System for the identification of yeasts

SUMMARY AND EXPLANATION

API Candida is a standardized system for the identification in 18-24 hours of yeasts, notably those most frequently encountered in clinical microbiology. The species which can be identified by the system are indicated in the Identification Table at the end of this package insert.

PRINCIPLE

The API Candida strip consists of 10 tubes containing dehydrated substrates, which enable the performance of 12 identification tests (sugar acidification or enzymatic reactions).

The reactions produced during incubation are revealed by spontaneous color changes.

The reactions are read visually according to the Reading Table and identification is obtained by consulting the list of profiles in this package insert or using the identification software.

CONTENT OF THE KIT (Kit for 10 tests)

- 10 API Candida strips
- 10 ampules of API NaCl 0.85 % Medium (2 ml)
- 10 incubation boxes
- 10 result sheets
- 1 package insert

COMPOSITION

Strip

The composition of the API Candida strip is given in the Reading Table of this package insert.

Medium

API NaCl 0.85 % Medium 2 ml	Sodium chloride Demineralized water	8.5 g 1000 ml
---	--	------------------

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents / Instrumentation


- McFarland Standard (ref. 70 900), point 3 on the scale
- Mineral oil (ref. 70 100)
- **apiweb**™ identification software (Ref. 40 011), (consult bioMérieux)

Material

- Pipettes or PSipettes
- Ampule protector
- Ampule rack
- Swabs
- General microbiology laboratory equipment


WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- All specimens, yeast cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling yeasts should be observed throughout this procedure. Refer to "CLSI/NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline - Current revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiry date.
- Before use, check that the packaging and components are intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.
- Allow media to come to room temperature before use.
- Open ampules carefully as follows :

* *Model 1 :*

- Cover the flattened part of the cap with the upper part of the thumb.
- Apply thumb pressure to the base of the flattened part of the cap to snap off the top of the ampule.

* *Model 2 :*

- Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.
- Take the ampule out of the ampule protector and put the protector aside for subsequent use.
- Carefully remove the cap.

- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed.

STORAGE CONDITIONS

Strips

The strips should be stored at 2-8°C until the expiry date indicated on the packaging.

Media

The media may be stored at 2-30°C until the expiry date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

API Candida is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Selection of colonies

Exam the strain to be studied under a microscope to check that it is indeed a yeast.

The colonies may be isolated on the following media before using the API Candida strip :

- Sabouraud 2 agar (with or without antibiotic) ;
- Albicans ID 2 agar ;
- Blood agar ;
- If other media are used to isolate the colonies, perform a subculture on one of the media mentioned above.

Preparation of the strip

- Prepare the incubation box (tray and lid) and distribute about 5 ml of water [demineralized or distilled water, or any water without additives or chemicals which may release gases (e.g. Cl₂, CO₂, etc.)] into the honey-combed wells of the tray to create a humid atmosphere.
- Record the strain reference on the elongated flap of the tray.
- Remove the strip from its individual packaging.
- Place the strip in the incubation box.
- Discard the desiccant sachet.

Preparation of the inoculum

- Open an ampule of API NaCl 0.85 % Medium (2 ml) as indicated in the paragraph "Warnings and Precautions".
- Using a pipette or swab, take one or several well-isolated, identical colonies and prepare a suspension with a turbidity equivalent to 3 McFarland : compare with a turbidity control or measure using a densitometer. It is recommended to use young cultures (18-24 hours old).
- Homogenize the yeast suspension. This suspension must be used immediately after preparation.

Inoculation of the strip

- Distribute the prepared yeast suspension into the tubes only, avoiding the formation of bubbles (tilt the incubation box slightly forward and position the pipette or PSIpette on the edge of the cupule).
- Cover the first 5 tests (GLU to RAF) and the last test (URE) with mineral oil (underlined tests) immediately after inoculating the strip.

NOTE : The quality of the filling is very important : tubes which are insufficiently or excessively full may cause false positive or false negative results.

- Close the incubation box.
- Incubate for 18-24 hours at 36°C ± 2°C in **aerobic conditions**.

READING AND INTERPRETATION

Reading the strip

After the 18-24-hours of incubation :


- Read the reactions by referring to the reading table in this package insert and record them as + or – on the result sheet.

Warning : Tubes 8 and 9 are bi-functional and enable 2 reactions to be performed in the same tube :

- Tube 8 : βXYL (test no. 8) / βNAG (test no. 11).
- Tube 9 : βGUR (test no. 9) / βGAL (test no. 12).

Interpretation

- Code the reactions obtained into a **numerical profile** : on the result sheet, the tests are separated into groups of three and a number 1, 2 or 4 is assigned to each one. By adding together the numbers corresponding to positive reactions within each group, a 4-digit numerical profile is obtained.
- Identification :
This is performed using the database (V 2.1)
* with the numerical profile :
- Look up the profile in the list in this package insert.
* with the **apiweb™** identification software :
- Enter the 4-digit numerical profile manually via the keyboard.
- If there is low discrimination between several species, complementary tests are indicated (see tables at the end of this package insert) to separate them. The results of these tests are taken from literature (or the ID 32 C database).

QUALITY CONTROL

The strips and media are systematically quality controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain 1. **Candida kefyr ATCC® 4135** or else one of the following strains :

- | | | | | |
|----|------------------------------|------|---------|--|
| 2. | <i>Trichosporon mucoides</i> | ATCC | 201382* | ATCC : American Type Culture Collection, |
| 3. | <i>Candida glabrata</i> | ATCC | 2001 | 10801 University Boulevard, Manassas, |
| | | | | VA 20110-2209, USA. |

	<u>GLU</u>	<u>GAL</u>	<u>SAC</u>	<u>TRE</u>	<u>RAF</u>	β MAL	α AMY	β XYL	β GUR	<u>URE</u>	β NAG	β GAL
1.	+	+	+	-	+	-	-	+	-	-	-	+
2.	+	+	+	-	V	+	+	-	+	+	+	-
3.	+	-	-	+	-	-	-	-	-	-	-	-

Profiles obtained after culture of the strains on Sabouraud agar.

* *Trichosporon mucoides* identified as *Trichosporon* spp 1 with API Candida.

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

- The API Candida system is designed uniquely for the identification of the species included in the database (see Identification Table at the end of this package insert) i.e. those belonging to the genera *Candida*, *Cryptococcus*, *Geotrichum*, *Saccharomyces* and *Trichosporon*. It cannot be used to identify any other microorganisms or to exclude their presence.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

PERFORMANCE

- 646 collection strains and strains of various origins belonging to species included in the database were tested:
 - 97.99 % of the strains were correctly identified (with or without supplementary tests).
 - 0.46 % of the strains were not identified.
 - 1.55 % of the strains were misidentified.

WASTE DISPOSAL

Unused reagents may be considered as non hazardous waste and disposed of accordingly.

Dispose of all used reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

READING TABLE

TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS/ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
1) <u>GLU</u>	D-glucose	1.4	Acidification (GLUcose)	violet grey-violet	yellow green / grey
2) <u>GAL</u>	D-galactose	1.4	Acidification (GALactose)		
3) <u>SAC</u>	D-saccharose	1.4	Acidification (SACcharose)		
4) <u>TRE</u>	D-trehalose	1.4	Acidification (TREhalose)		
5) <u>RAF</u>	D-raffinose	1.4	Acidification (RAFfinose)		
6) <u>βMAL</u>	4-nitrophenyl-βD-maltopyranoside	0.08	β-MALtosidase	colorless	pale yellow-bright yellow
7) <u>αAMY</u>	2-chloro-4-nitrophenyl-αD-maltotrioside	0.168	α-AMYlase	colorless	pale yellow-bright yellow
8) <u>βXYL</u>	4-nitrophenyl-βD-xylopyranoside	0.095	β-XYLosidase	colorless-very pale yellow / blue / green **	pale yellow-bright yellow
9) <u>βGUR</u>	4-nitrophenyl-βD-glucuronide	0.063	β-GIUcuRonidase	colorless / blue / green	pale yellow-bright yellow
10) <u>URE</u>	urea	1.68	UREase	yellow-pale orange	red
11) <u>βNAG</u> (in tube no. 8) *	5-bromo-4-chloro-3-indoxyl-N-acetyl-βD-glucosaminide	0.09	N-Acetyl-β-Glucosaminidase	colorless / yellow	blue / green **
12) <u>βGAL</u> (in tube no. 9) *	5-bromo-4-chloro-3-indolyl-βD-galactopyranoside	0.0815	β-GALactosidase	colorless / yellow	blue / green

* Tubes 8 and 9 are bi-functional : tube 8 : βXYL (test no. 8) / βNAG (test no. 11)
tube 9 : βGUR (test no. 9) / βGAL (test no. 12)

** Any trace of green in cupule 8 = βXYL (-) βNAG (+)

- The quantities indicated may be adjusted depending on the titer of the raw materials used.

PROCEDURE	p. I
LIST OF NUMERICAL PROFILES	p. II
TABLES	p. III
IDENTIFICATION TABLE	p. V
LITERATURE REFERENCES	p. VII
INDEX OF SYMBOLS	p. VIII

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

bioMérieux® SA
au capital de 12 029 370 €
673 620 399 RCS LYON
69280 Marcy-l'Etoile / France
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

bioMérieux, Inc
Box 15969,
Durham, NC 27704-0969 / USA
Tel. (1) 919 620 20 00
Fax (1) 919 620 22 11

Printed in France

bioMérieux, the blue logo, API and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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The following morphological characteristics should be noted /

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1. BARNETT J.A., PAYNE R.W., YARROW D.
Yeasts : Characteristics and Identification.
(1990) Cambridge University Press, London.
2. BERNAL S., MARTIN MAZUELOS E., CHAVEZ M.,
CORONILLA J., VALVERDE A.
Evaluation of the new API Candida system for identification
of the most clinically important yeast species.
(1998) Diagn. Microbiol. Infect. Dis., 32, 3, 217-221.
3. DURUSSEL C., BILLE J.
API Candida, a New Simplified 12 Tests Rapid Identification
System for Yeasts.
(1996) 96th ASM General Meeting, New Orleans, Louisiana,
F-78.
4. FRICKER-HIDALGO H., VANDAPEL O., DUCHESNE M.A.,
MAZOYER M.A., MONGET D., LARDY B., LEBEAU B.,
FRENEY J., AMBROISE-THOMAS P., GRILLOT R.
Comparison of the New API Candida System to the ID 32C
System for the Identification of Clinically Important Yeast
Species.
(1996) J. Clin. Microbiol., 34, 1846-1848.
5. KREGER VAN RIJ N.J.W.
The Yeasts : A Taxonomic Study.
(1984) Elsevier, Amsterdam.
6. LARONE D.H.
Medically Important Fungi. A Guide to Identification.
Third Edition.
(1995) A.S.M., Washington, D.C.
7. MCGINNIS M.R. and al.
Taxonomic and Nomenclatural Evaluation of the genera
Candida and *Torulopsis*.
(1984) J. Clin. Microbiol., 20, 813-814.
8. MONGET D., DUCHESNE M.A., CANIAUX I.
api Candida, A New Identification System for Yeasts.
(1995) 7th E.C.C.M.I.D., Vienna, 26-30 March 1995.
9. WARREN N.G., HAZEN K.C.
Candida, *Cryptococcus*, and Other Yeasts of Medical
Importance in Manual of Clinical Microbiology.
Seventh Edition.
(1995) A.S.M., Washington, D.C., 95, 1184-1199.

**TABLE DES SYMBOLES / INDEX OF SYMBOLS / SIMBOLE / CUADRO DE SIMBOLOS /
TABELLA DEI SIMBOLI / QUADRO DOS SÍMBOLOS / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩΝ / SYMBOLER /
SYMBOLFORTEGNELSE / TABELA SYMBOLI**

Symbole / Symbol Símbolo / Simbolo Σύμβολο	Signification / Meaning / Bedeutung Significado / Significato / Επεξήγηση Betydelse / Betydning / Znaczenie
	Référence du catalogue Catalogue number (GB) / Catalog number (US) Bestellnummer / Número de catálogo / Numero di catalogo Referência de catálogo / Αριθμός καταλόγου Katalognummer / Katalognummer / Numer katalogowy
	Dispositif médical de diagnostic in vitro In Vitro Diagnostic Medical Device / In Vitro Diagnostikum Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro Dispositivo médico para diagnóstico in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Medicintekniska produkter för in vitro diagnostik Medicinsk udstyr til in vitro-diagnostik Wyrób do diagnostyki In Vitro
	Fabricant / Manufacturer / Hersteller / Fabricante Fabbicante / Κατασκευαστής / Tillverkare / Producent
	Limites de température / Temperature limitation Temperaturbegrenzung / Limite de temperatura Limiti di temperatura / Limites de temperatura Περιορισμοί θερμοκρασίας / Temperaturbegränsning Temperaturbegrænsning / Przestrzegać zakresu temperatury
	Utiliser jusque / Use by / Verwendbar bis Fecha de caducidad / Utilizzare entro / Prazo de validade Ημερομηνία λήξης / Använd före / Holdbar til / Użyć przed
	Code du lot / Batch code / Chargenbezeichnung Código de lote / Codice del lotto / Código do lote Αριθμός Παρτίδας / Lot nummer / Lotnummer / Kod partii
	Consulter les instructions d'utilisation Consult Instructions for Use / Gebrauchsanweisung beachten Consulte las instrucciones de uso Consultare le istruzioni per l'uso Consulte as instruções de utilização Συμβουλευτείτε τις οδηγίες χρήσης Se handhavandebeskrivningen / Se brugsanvisning Sprawdź w instrukcji obsługi
	Contenu suffisant pour "n" tests Contains sufficient for <n> tests Inhalt ausreichend für <n> Prüfungen Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi Conteúdo suficiente para "n" ensaios Περιεχόμενο επαρκές για «n» εξετάσεις Räcker till "n" antal tester Indeholder tilstrækkeligt til "n" test Wystarczy na wykonanie <n> testów